



an-AlayVet

Cross-Linked HA for Lid Augmentation

- CIS Technology. This “core-in-shell” structure forms a stable core of highly cross-linked hyaluronic acid, surrounded by long chains of non-cross-linked hyaluronic acid. This technology makes it different from other products made of cross-linked hyaluronic acid.
- Delayed degradation of hyaluronic acid by hyaluronidase
- Prolonged effects of the treatment
- Skin filler for effective treatment of entropion
- Free of animal components
- Excellent biocompatibility (clinically proven)
- No allergy risk





Photo Credit: Dr. Amy Knollinger



Directions:

- Disinfect injection site
- Apply local anesthetic as needed
- Use **an-AlayVet** at room temperature
- Insert the included 27G needle into the medium or deep dermal level
- Inject slowly while simultaneously backing out the needle from the skin.

Specifications

- Sodium Hyaluronate 24mg / 1ml
- Viscosity 30.000-50.000 mPa.s
- pH 6.8 - 7.3
- Molecular Weight 0.8 - 1.2 x 10⁶ Dalton
- Packaging Box contains 2 syringes, with 1ml in each



an-vision GmbH
T: +49 (0)3302-209 5550
info@an-vision.eu

www.an-vision.com

an-vision Inc.
T: 801-561-5040
info@an-vision.com

Day 0
Pre-Injection



Day 0
Post-Injection



Day 21
Follow-Up



Day 60
Follow-Up



Photos Courtesy of Dr. Jose Sanabria