

## Bandage Lens Sizing Chart

Bandage Lens Sizes for Cats			
	Diameter	Base Curve	
C1	16.0 mm	9.2 mm	Cats normal - slightly larger
C2	14.0 mm	8.3 mm	Kitten
C3	15.5 mm	8.5 mm	Young Cats
C4	16.0 mm	8.8 mm	Cats normal - slightly smaller
C5	18.0 mm	9.4 mm	Persian Cats - large eyes
Bandage Lens Sizes for Dogs			
	Diameter	Base Curve	
D1	18.0 mm	9.6 mm	Medium sized dogs - somewhat smaller; most dog breeds such as collies, huskies, corgis, some boxers
D2	18.0 mm	9.8 mm	Medium sized dogs - somewhat larger; most dog breeds such labradors, retrievers, boxers
D3	20.0 mm	10.4 mm	Dogs - large eyes; some bulldogs; some boxers
D4	15.5 mm	8.5 mm	Young dogs
D5	22.0 mm	11.4 mm	Dogs - very large eyes
D6	22.0 mm	11.8 mm	Dogs - very large eyes only
D7	20.0 mm	11.2 mm	Dogs - in between sizes
D8	20.0 mm	10.8 mm	Dogs - in between sizes
D9	18.0 mm	10.0 mm	Dogs - in between sizes
Bandage Lens Sizes for Horses			
	Diameter	Base Curve	
H1	34.0 mm	18.0 mm	Horses, Normal Eye Size
H2	32.0 mm	17.0 mm	Small Horses, Ponies
H3	26.0 mm	15.0 mm	Young Ponies
H4	30.0 mm	16.0 mm	Young Horses
H5	38.0 mm	19.5mm	Large Horses
H6	36.0 mm	20.0 mm	Horses, Larger Eyes