


an-vision

for animal vision

an-Lens UPGRADE

Take Advantage of the MD8 Design



The **an-lens MD6** and **MD8** have the:


- **same optic size**
- **same haptic design**
- **same optical design**

The **MD8** offers the added advantage of:

- **the additional 360° barrier to prevent PCO**

Use of the MD8

- To take advantage of the posterior barrier, simply load the lens inside the cartridge using the clearly marked notches on the haptics
- If the lens is implanted the other way, it remains identical to the MD6



We offer the advanced MD8 for the **same price** as the MD6!

an-vision Inc.
T: 801-561-5040
info@an-vision.com

www.an-vision.com

an-vision GmbH
T: +49 (0)3302-209 5550
info@an-vision.eu